


Bioline™

HEPATITIS B SERIES

HEPATITIS B VIRUS TEST

HBsAg, HBsAg WB

Bioline™ Hepatitis tests are intended for professional use as an aid in the diagnosis of hepatitis B. Highly sensitive, specific immunochromatographic assays for detection of HBsAg.




	HBsAg	HBsAg WB
SPECIMEN	Serum, Plasma	Serum, Plasma, Whole blood
SENSITIVITY	100 %	100 %
SPECIFICITY	100 %	100 %

MATERIALS PROVIDED


- Test device

SIMPLE PROCEDURE


- Add Specimen**
Dispense 100 µl of specimen into the specimen well.



Wait 20 mins.



RESULTS INTERPRETATION



The presence of any test line, no matter how faint, the result is considered positive.

ORDERING INFORMATION

PRODUCT	CAT. NO.	TYPE	SPECIMEN	PACK SIZE
HBsAg	01FK10	Device	Serum/Plasma	30T/Kit
HBsAg WB	01FK10W	Device	Serum/Plasma/Whole blood	30T/Kit